

PAXIP1 Antibody



PACO11085

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAXIP1

Uniprot

Q6ZW49

Synonyms:

PAX interacting (with transcription-activation domain) protein 1; PAXIP1; CAGF28; CAGF29; FLJ41049; PACIP1; PAXIP1L; PTIP; TNRC2 ;

Immunogen:

Human PAXIP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A