

# PBXIP1 Antibody



PACO11087

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PBXIP1

**Uniprot**

Q96AQ6

**Synonyms:**

pre-B-cell leukemia homeobox interacting protein 1; PBXIP1; HPIP ;

**Immunogen:**

Human PBXIP1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A