

PCBP4 Antibody



PACO11091

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PCBP4

Uniprot

P57723

Synonyms:

42kd;

Immunogen:

Human PCBP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A