

PCDH11X Antibody



PACO11093

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PCDH11X

Uniprot

Q9BZA7

Synonyms:

protocadherin 11 X-linked; PCDH11X; PCDH-X; PCDH11; PCDHX; PCDHY ;

Immunogen:

Human PCDH11X.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A