

PCDH11Y Antibody



PACO11094

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCDH11Y

Uniprot

Q9BZA8

Synonyms:

protocadherin 11 Y-linked; PCDH11Y; PCDH-PC; PCDH11X; PCDH22; PCDHY ;

Immunogen:

Human PCDH11Y.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A