

# PCDHA2 Antibody



PACO11095

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCDHA2

**Uniprot**

Q9Y5H9

**Synonyms:**

protocadherin alpha 2; PCDHA2; PCDH-ALPHA2 ;

**Immunogen:**

Human PCDHA2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A