

# PCDHA3 Antibody



PACO11096

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCDHA3

**Uniprot**

Q9Y5H8

**Synonyms:**

protocadherin alpha 3; PCDHA3; MGC141669; PCDH-ALPHA3 ;

**Immunogen:**

Human PCDHA3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A