

PCDHA5 Antibody



PACO11097

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCDHA5

Uniprot

Q9Y5H7

Synonyms:

protocadherin alpha 5; PCDHA5; CNR6; CNRN6; CNRS6; CRNR6; PCDH-ALPHA5 ;

Immunogen:

Human PCDHA5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A