

PCDHA6 Antibody



PACO11098

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCDHA6

Uniprot

Q9UN73

Synonyms:

protocadherin alpha 6; PCDHA6; CNR2; CNRN2; CNRS2; CRNR2; PCDH-ALPHA6 ;

Immunogen:

Human PCDHA6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A