

PCDHB12 Antibody



PACO11099

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCDHB12

Uniprot

Q9Y5F1

Synonyms:

protocadherin beta 12; PCDHB12; PCDH-BETA12 ;

Immunogen:

Human PCDHB12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A