

PCDHB5 Antibody



PACO11101

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCDHB5

Uniprot

Q9Y5E4

Synonyms:

protocadherin beta 5; PCDHB5; DKFZp586B0217; PCDH-BETA5 ;

Immunogen:

Human PCDHB5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A