

PCDHGB3 Antibody



PACO11102

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCDHGB3

Uniprot

Q9Y5G1

Synonyms:

protocadherin gamma subfamily B, 3; PCDHGB3; PCDH-GAMMA-B3 ;

Immunogen:

Human PCDHGB3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A