

PCK2 Antibody



PACO11106

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCK2

Uniprot

Q16822

Synonyms:

phosphoenolpyruvate carboxykinase 2 (mitochondrial); PCK2; PEPCK; PEPCK-M; PEPCK2 ;

Immunogen:

Human PCK2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A