

# PCMT1 Antibody



PACO11107

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCMT1

**Uniprot**

P22061

**Synonyms:**

protein-L-isoaspartate (D-aspartate) O-methyltransferase;

**Immunogen:**

Human PCMT1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A