

# PCNA Antibody



PACO11108

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat, Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCNA

**Uniprot**

P12004

**Synonyms:**

proliferating cell nuclear antigen; PCNA; MGC8367 ;

**Immunogen:**

Human PCNA.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A