

PCNP Antibody



PACO11109

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCNP

Uniprot

Q8WW12

Synonyms:

PEST proteolytic signal containing nuclear protein; PCNP; DKFZp781I24156 ;

Immunogen:

Human PCNP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A