

PCNXL2 Antibody



PACO11110

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PCNX2

Uniprot

A6NKB5

Synonyms:

pecanex-like 2 (Drosophila); PCNXL2; FLJ11383; KIAA0435 ;

Immunogen:

Human PCNXL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A