

PCSK4 Antibody



PACO11113

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCSK4

Uniprot

Q6UW60

Synonyms:

proprotein convertase subtilisin/kexin type 4; PCSK4; DKFZp434B217; MGC34749; PC4; SPC5 ;

Immunogen:

Human PCSK4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A