

PCYT2 Antibody



PACO11122

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCYT2

Uniprot

Q99447

Synonyms:

phosphate cytidyltransferase 2, ethanolamine; PCYT2; ET ;

Immunogen:

Human PCYT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A