

PDCD7 Antibody



PACO11131

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PDCD7

Uniprot

Q8N8D1

Synonyms:

programmed cell death 7; PDCD7; ES18; HES18; MGC22015 ;

Immunogen:

Human PDCD7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A