

PDCL Antibody



PACO11132

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PDCL

Uniprot

Q13371

Synonyms:

phosducin-like; PDCL; DKFZp564M1863; PhLP ;

Immunogen:

Human PDCL

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A