

PDGFD Antibody



PACO11134

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PDGFD

Uniprot

Q9GZP0

Synonyms:

platelet derived growth factor D; PDGFD; IEGF; MGC26867; MSTP036; SCDGF-B; SCDGFB ;

Immunogen:

Human PDGFD.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A