

PDIA3 Antibody



PACO11139

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

PDIA3

Uniprot

P30101

Synonyms:

protein disulfide isomerase family A, member 3; PDIA3; ER60; ERp57; ERp60; ERp61; GRP57; GRP58; HsT17083; P58; PI-PLC ;

Immunogen:

Human PDIA3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A