

PDIA5 Antibody



PACO11141

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PDIA5

Uniprot

Q14554

Synonyms:

protein disulfide isomerase family A, member 5; PDIA5; FLJ30401; PDIR ;

Immunogen:

Human PDIA5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A