

PDIA6 Antibody



PACO11142

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PDIA6

Uniprot

Q15084

Synonyms:

protein disulfide isomerase family A, member 6; PDIA6; ERP5; P5; TXNDC7 ;

Immunogen:

Human PDIA6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A