

PDLIM3 Antibody



PACO11146

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PDLIM3

Uniprot

Q53GG5

Synonyms:

PDZ and LIM domain 3; PDLIM3; ALP; DKFZp686L0362 ;

Immunogen:

Human PDLIM3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A