

PDLIM7 Antibody



PACO11148

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PDLIM7

Uniprot

Q9NR12

Synonyms:

PDZ and LIM domain 7 (enigma); PDLIM7; LMP1 ;

Immunogen:

Human PDLIM7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A