

PDRG1 Antibody



PACO11150

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

PDRG1

Uniprot

Q9NUG6

Synonyms:

p53 and DNA damage regulated 1; PDRG1; C20orf126; PDRG ;

Immunogen:

Human PDRG1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A