

PDXDC1 Antibody



PACO11154

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PDXDC1

Uniprot

Q6P996

Synonyms:

pyridoxal-dependent decarboxylase domain containing 1; PDXDC1; KIAA0251; LP8165 ;

Immunogen:

Human PDXDC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A