

ECI2 Antibody



PACO11161

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ECI2

Uniprot

O75521

Synonyms:

peroxisomal D3, D2-enoyl-CoA isomerase; PECl; ACBD2; DRS1; HCA88; KIAA0536; dJ1013A10.3 ;

Immunogen:

Human PECl.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A