

PES1 Antibody



PACO11171

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PES1

Uniprot

O00541

Synonyms:

pescadillo homolog 1, containing BRCT domain (zebrafish); PES1; PES ;

Immunogen:

Human PES1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A