

PFDN2 Antibody



PACO11180

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PFDN2

Uniprot

Q9UHV9

Synonyms:

prefoldin subunit 2; PFDN2; PFD2 ;

Immunogen:

Human PFDN2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A