

# PFKFB2 Antibody



PACO11183

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PFKFB2

**Uniprot**

O60825

**Synonyms:**

6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 2; PFKFB2; DKFZp781D2217; MGC138308; MGC138310; PFK-2/FBPase-2 ;

**Immunogen:**

Human PFKFB2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A