

# PFKP Antibody



PACO11187

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PFKP

**Uniprot**

Q01813

**Synonyms:**

phosphofructokinase, platelet; PFKP; FLJ40226; PFK-C; PFKF ;

**Immunogen:**

Human PFKP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A