

PGBD3 Antibody



PACO11192

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PGBD3

Uniprot

Q8N328

Synonyms:

piggyBac transposable element derived 3; PGBD3; FLJ90201 ;

Immunogen:

Human PGBD3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A