

PACO11195

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PGF

**Uniprot**

P49763

**Synonyms:**

placental growth factor; PGF; D12S1900; PGFL; PLGF; PIGF-2; SHGC-10760 ;

**Immunogen:**

Human PGF.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A