

PGK2 Antibody



PACO11197

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PGK2

Uniprot

P07205

Synonyms:

phosphoglycerate kinase 2; PGK2; PGK-2; PGKB; PGKPS; dJ417L20.2 ;

Immunogen:

Human PGK2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A