

PHACTR4 Antibody



PACO11205

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PHACTR4

Uniprot

Q8IZ21

Synonyms:

phosphatase and actin regulator 4; PHACTR4; DKFZp686L07205; FLJ13171; MGC20618; MGC34186 ;

Immunogen:

Human PHACTR4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A