

# PHB2 Antibody



PACO11208

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PHB2

**Uniprot**

Q99623

**Synonyms:**

prohibitin 2; PHB2; BAP; BCAP37; Bap37; MGC117268; PNAS-141; REA; p22 ; prohibitin-2

**Immunogen:**

Human PHB2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A