

# PHIP Antibody



PACO11215

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PHIP

**Uniprot**

Q8WWQ0

**Synonyms:**

pleckstrin homology domain interacting protein; PHIP; FLJ20705; FLJ45918; MGC90216; WDR11; ndrp ;

**Immunogen:**

Human PHIP.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A