

PHLDA1 Antibody



PACO11219

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PHLDA1

Uniprot

Q8WV24

Synonyms:

pleckstrin homology-like domain, family A, member 1; PHLDA1; DT1P1B11; MGC131738; PHRIP; TDAG51 ;

Immunogen:

Human PHLDA1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A