

PHPT1 Antibody



PACO11221

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PHPT1

Uniprot

Q9NRX4

Synonyms:

phosphohistidine phosphatase 1; PHPT1; CGI-202; DKFZp564M173; HSPC141; PHP14; bA216L13.10 ;

Immunogen:

Human PHPT1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A