

PHYH Antibody



PACO11222

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PHYH

Uniprot

O14832

Synonyms:

phytanoyl-CoA 2-hydroxylase; PHYH; LN1; LNAP1; PAHX; PHYH1; RD ;

Immunogen:

Human PHYH.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A