

# PHYHIPL Antibody



PACO11225

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PHYHIPL

**Uniprot**

Q96FC7

**Synonyms:**

phytanoyl-CoA 2-hydroxylase interacting protein-like; PHYHIPL; KIAA1796 ;

**Immunogen:**

Human PHYHIPL.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A