

PIGB Antibody



PACO11237

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PIGB

Uniprot

Q92521

Synonyms:

phosphatidylinositol glycan anchor biosynthesis, class B; PIGB; MGC21236 ;

Immunogen:

Human PIGB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A