

PIGK Antibody



PACO11238

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PIGK

Uniprot

Q92643

Synonyms:

phosphatidylinositol glycan anchor biosynthesis, class K; PIGK; GPI8; MGC22559 ;

Immunogen:

Human PIGK.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A