

# PIGM Antibody



PACO11239

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PIGM

**Uniprot**

Q9H3S5

**Synonyms:**

phosphatidylinositol glycan anchor biosynthesis, class M; PIGM; GPI-MT-I; MGC29896 ;

**Immunogen:**

Human PIGM.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A