

PIGO Antibody



PACO11240

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PIGO

Uniprot

Q8TEQ8

Synonyms:

phosphatidylinositol glycan anchor biosynthesis, class O; PIGO; DKFZp434M222; FLJ00135; MGC20536; MGC3079; RP11-182N22.4 ;

Immunogen:

Human PIGO.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A