

PIGZ Antibody



PACO11242

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PIGZ

Uniprot

Q86VD9

Synonyms:

phosphatidylinositol glycan anchor biosynthesis, class Z; PIGZ; FLJ12768; GPI-MT-IV; MGC52163; SMP3 ;

Immunogen:

Human PIGZ.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A