

PIP4K2A Antibody



PACO11253

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PIP4K2A

Uniprot

P48426

Synonyms:

phosphatidylinositol-5-phosphate 4-kinase, type II, alpha; PIP4K2A; FLJ13267; PI5P4KA; PIP5K2A; PIP5KII-alpha; PIP5KIIA; PIPK ;

Immunogen:

Human PIP4K2A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A